

Bed Grab Rail



Designed specifically to fit between the mattress and base of a divan bed, this high quality grab rail aids the elevation of the body from horizontal to a seated position and also assists in transfer to a wheelchair, commode or standing position.

A nylon strap and buckle enables the rail to be fixed securely to the bed base.

The metal frame is coated with a magnolia paint containing an antibacterial agent, which helps prevent the transfer of infection.

Free-phone for any queries or a delivery date on 0800 194 007

**1 YEAR
GUARANTEE**

**7 DAY
MONEYBACK
GUARANTEE**

**FREE
UK MAINLAND
DELIVERY**